IT SUNEKOS PERFORMA

Performa

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SUNEKOS PERFORMA

Performa

PATIENT INFORMATION LEAFLET**DESCRIPTION**

Sunekos® Performa is an implantable and resorbable device based on non-crosslinked (linear) hyaluronic acid of non-animal origin and a specific mixture of six amino acids.

Sunekos® Performa is packaged in single-use doses.

CATEGORIES OF USERS AND CONSUMERS (Who can use it and who is it intended for?)

Sunekos® Performa should only be administered by healthcare professionals with appropriate training and experience.

Treatment with Sunekos® Performa is intended for adults with no gender differences and aged 18 years and over.

EXPECTED RESULTS (What does Sunekos® Performa do?)

Sunekos® Performa promotes the correction of imperfections and skin depressions of the face and neck of mild to moderate severity.

In particular the product:

- exert a filling effect on anatomically depressed areas of the skin, such as wrinkles and folds, promoting the normalization of imperfections and,
- promote the restoration of hydration and the physiological elasticity of the skin tissue.

CONTRAINDICATIONS (Are there any reasons why you should not undergo treatment with Sunekos® Performa?)

You should not undergo treatment with Sunekos® Performa if:

- You are allergic or hypersensitive to the components (sodium hyaluronate and amino acids (glycine, proline, leucine, lysine, alanine, valine));
- You are under 18 years old;
- You are pregnant or breastfeeding;
- You are suffering from inflammation and/or active infections;
- You are undergoing anticoagulant therapy;
- You are suffering from autoimmune diseases, neoplastic diseases (tumors) and/or undergoing immunosuppressant and chemotherapy treatments.

PRECAUTIONS FOR USE

Before starting the treatment, the patient should inform the healthcare professional of any past or current conditions, chronic diseases, allergies and current therapies, as well as any previous treatments with other permanent or temporary implants so that the healthcare professional can take an accurate medical history. Before the injection, the healthcare professional carrying out the treatment must provide the patient with the informed consent and this information leaflet, containing all the information useful to the patient, such as indications, warnings, contraindications, and potential known side effects.

After the treatment, the consumer/patient must receive the Patient Implant Card (PIC) from the healthcare professional who performed the injection.

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The healthcare professional is required to fill in the PIC with all the information required by the manufacturer and useful to identify the device (batch and UDI) and the person who carried out the treatment.
The patient is required to keep the PIC carefully.

WARNINGS

The product is single-use/single-patient and should only be used in a single session, reuse of the product for separate applications is explicitly prohibited and may cause infection for the consumer/patient.

It is recommended that the consumer/patient does not expose the treated area to heat sources until the swelling has been reabsorbed.

The product should not be injected intramuscularly or intravascularly, the introduction of the latter into the vascular system may lead to embolism, vessel occlusion, ischemia or infarction.

Rare but serious adverse events associated with intravascular injection of soft tissue fillers containing hyaluronic acid into the face are known in the literature, including temporary or permanent visual disturbances, blindness, cerebral ischemia or cerebral haemorrhage, resulting in stroke, skin necrosis and damage to underlying facial structures.

The injection should be stopped immediately if the consumer/patient has any of the following symptoms, including changes in vision, signs of stroke, pale skin, or unusual pain during or shortly after the procedure. If an intravascular injection occurs, the patient should receive immediate medical attention and possibly an evaluation by an appropriate specialist physician.

Please note that the aforementioned adverse events have not been reported in relation to the administration of Sunekos® Performa.

POSSIBLE SIDE EFFECTS

Sunekos® Performa is characterized by high tolerability. However, as with all hyaluronic acid injectable products, local reactions may occur after injection, which may manifest as pain, sensations of heat and/or itching, swelling, inflammatory reactions (redness), bleeding and bruising. Reactions usually resolve within a few days.

The device does not contain carcinogenic, mutagenic or reprotoxic substances.

In the event of undesirable adverse events not reported in this patient leaflet or if the reported effects persist beyond the indicated duration, please contact the person who carried out the treatment and/or the manufacturer Professional Dietetics S.p.A. through the following contacts:

Phone +39.02.74.40.20

info@pd-international.com

Serious incidents that should be reported to the Therapeutic Goods Administration (TGA) Australia

<https://www.tga.gov.au>

This patient leaflet is in compliance with Annexes I and IV of Implementing Regulation (EU) 2022/2346.

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Soggetta alla Direzione e Coordinamento di Iera S.r.l.